

Anti-CFAP91 antibody (310-390) (STJ190542)

STJ190542

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cilia-And Flagella-Associated Protein 91 (310-390) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89876
Gene Symbol	CFAP91
Uniprot ID	CFA91_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 310-390
Immunogen Region	310-390
Specificity	CFAP91 polyclonal antibody (Cilia-And Flagella-Associated Protein 91) binds to endogenous Cilia-And Flagella-Associated Protein 91 at the amino acid region 310-390.
Immunogen Sequence	