

Anti-RECQL4 antibody (1030-1110) (STJ190591)

STJ190591

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Atp-Dependent Dna Helicase Q4 (1030-1110) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9401
Gene Symbol	RECQL4
Uniprot ID	RECQ4_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1030-1110
Immunogen Region	1030-1110
Specificity	RECQL4 polyclonal antibody (Atp-Dependent Dna Helicase Q4) binds to endogenous Atp-Dependent Dna Helicase Q4 at the amino acid region 1030-1110.
Immunogen Sequence	