

Anti-PDE3B antibody (STJ190687)

STJ190687

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cgmp-Inhibited 3 5 Cyclic Phosphodiesterase B is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5140
Gene Symbol	PDE3B
Uniprot ID	PDE3B_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	PDE3B polyclonal antibody (Cgmp-Inhibited 3 NA-5 NA-Cyclic Phosphodiesterase B) binds to endogenous Cgmp-Inhibited 3 NA-5 NA-Cyclic Phosphodiesterase B.
Immunogen Sequence	