

Anti-PPFIA1 antibody (980-1060) (STJ190927)

STJ190927

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Liprin-Alpha-1 (980-1060) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8500
Gene Symbol	PPFIA1
Uniprot ID	LIPA1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 980-1060
Immunogen Region	980-1060
Specificity	PPFIA1 polyclonal antibody (Liprin-Alpha-1) binds to endogenous Liprin-Alpha-1 at the amino acid region 980-1060.
Immunogen Sequence	