

Anti-ACSL5 antibody (460-540) (STJ190933)

STJ190933

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Long-Chain-Fatty-Acid--CoA Ligase 5 (460-540) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51703
Gene Symbol	ACSL5
Uniprot ID	ACSL5_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 460-540
Immunogen Region	460-540
Specificity	ACSL5 polyclonal antibody (Long-Chain-Fatty-Acid--CoA Ligase 5) binds to endogenous Long-Chain-Fatty-Acid--CoA Ligase 5 at the amino acid region 460-540.
Immunogen Sequence	