

## Anti-LRP4 antibody (290-370) (STJ190937)

STJ190937

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Low-Density Lipoprotein Receptor-Related Protein 4 (290-370) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Formulation</b>	PBS, 50% Glycerol and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">4038</a>
<b>Gene Symbol</b>	<a href="#">LRP4</a>
<b>Uniprot ID</b>	<a href="#">LRP4_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from human protein at aa range 290-370
<b>Immunogen Region</b>	290-370
<b>Specificity</b>	LRP4 polyclonal antibody (Low-Density Lipoprotein Receptor-Related Protein 4) binds to endogenous Low-Density Lipoprotein Receptor-Related Protein 4 at the amino acid region 290-370.
<b>Immunogen Sequence</b>	