

Anti-LSM14A antibody (330-410) (STJ190938)

STJ190938

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Lsm14 Homolog A (330-410) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26065
Gene Symbol	LSM14A
Uniprot ID	LS14A_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 330-410
Immunogen Region	330-410
Specificity	LSM14A polyclonal antibody (Protein Lsm14 Homolog A) binds to endogenous Protein Lsm14 Homolog A at the amino acid region 330-410.
Immunogen Sequence	