

Anti-MYH8 antibody (760-840) (STJ191004)

STJ191004

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Myosin-8 (760-840) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4626
Gene Symbol	MYH8
Uniprot ID	MYH8_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 760-840
Immunogen Region	760-840
Specificity	MYH8 polyclonal antibody (Myosin-8) binds to endogenous Myosin-8 at the amino acid region 760-840.
Immunogen Sequence	