

Anti-MYO1F antibody (640-720) (STJ191009)

STJ191009

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Unconventional Myosin-I α (640-720) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4542
Gene Symbol	MYO1F
Uniprot ID	MYO1F_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 640-720
Immunogen Region	640-720
Specificity	MYO1F polyclonal antibody (Unconventional Myosin-I α) binds to endogenous Unconventional Myosin-I α at the amino acid region 640-720.
Immunogen Sequence	