

Anti-UNC5B antibody (350-430) (STJ191049)

STJ191049

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Netrin Receptor Unc5b (350-430) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	219699
Gene Symbol	UNC5B
Uniprot ID	UNC5B_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 350-430
Immunogen Region	350-430
Specificity	UNC5B polyclonal antibody (Netrin Receptor Unc5b) binds to endogenous Netrin Receptor Unc5b at the amino acid region 350-430.
Immunogen Sequence	