

Anti-NUP155 antibody (840-920) (STJ191096)

STJ191096

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nuclear Pore Complex Protein Nup155 (840-920) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9631
Gene Symbol	NUP155
Uniprot ID	NUP155_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 840-920
Immunogen Region	840-920
Specificity	NUP155 polyclonal antibody (Nuclear Pore Complex Protein Nup155) binds to endogenous Nuclear Pore Complex Protein Nup155 at the amino acid region 840-920.
Immunogen Sequence	