

Anti-PLEKHG1 antibody (780-860) (STJ191169)

STJ191169

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pleckstrin Homology Domain-Containing Family G Member 1 (780-860) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57480
Gene Symbol	PLEKHG1
Uniprot ID	PKHG1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 780-860
Immunogen Region	780-860
Specificity	PLEKHG1 polyclonal antibody (Pleckstrin Homology Domain-Containing Family G Member 1) binds to endogenous Pleckstrin Homology Domain-Containing Family G Member 1 at the amino acid region 780-860.
Immunogen Sequence	