

Anti-PCF11 antibody (490-570) (STJ191199)

STJ191199

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pre-Mrna Cleavage Complex 2 Protein Pcf11 (490-570) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51585
Gene Symbol	PCF11
Uniprot ID	PCF11_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 490-570
Immunogen Region	490-570
Specificity	PCF11 polyclonal antibody (Pre-Mrna Cleavage Complex 2 Protein Pcf11) binds to endogenous Pre-Mrna Cleavage Complex 2 Protein Pcf11 at the amino acid region 490-570.
Immunogen Sequence	