

Anti-PRPF40A antibody (450-530) (STJ191201)

STJ191201

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Pre-Mrna-Processing Factor 40 Homolog A (450-530) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55660
Gene Symbol	PRPF40A
Uniprot ID	PR40A_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 450-530
Immunogen Region	450-530
Specificity	PRPF40A polyclonal antibody (Pre-Mrna-Processing Factor 40 Homolog A) binds to endogenous Pre-Mrna-Processing Factor 40 Homolog A at the amino acid region 450-530.
Immunogen Sequence	