

Anti-PPM1E antibody (460-540) (STJ191215)

STJ191215

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Phosphatase 1e (460-540) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22843
Gene Symbol	PPM1E
Uniprot ID	PPM1E_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 460-540
Immunogen Region	460-540
Specificity	PPM1E polyclonal antibody (Protein Phosphatase 1e) binds to endogenous Protein Phosphatase 1e at the amino acid region 460-540.
Immunogen Sequence	