

Anti-RAB11FIP5 antibody (570-650) (STJ191237)

STJ191237

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Rab11 Family-Interacting Protein 5 (570-650) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26056
Gene Symbol	RAB11FIP5
Uniprot ID	RFIP5_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 570-650
Immunogen Region	570-650
Specificity	RAB11FIP5 polyclonal antibody (Rab11 Family-Interacting Protein 5) binds to endogenous Rab11 Family-Interacting Protein 5 at the amino acid region 570-650.
Immunogen Sequence	