

Anti-RFTN2 antibody (370-450) (STJ191244)

STJ191244

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Raftlin-2 (370-450) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	130132
Gene Symbol	RFTN2
Uniprot ID	RFTN2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 370-450
Immunogen Region	370-450
Specificity	RFTN2 polyclonal antibody (Raftlin-2) binds to endogenous Raftlin-2 at the amino acid region 370-450.
Immunogen Sequence	