

Anti-RGPD5/RGPD6 antibody (1560-1640) (STJ191249)

STJ191249

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ranbp2-Like And Grip Domain-Containing Protein 5/6 (1560-1640) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	729540/84220
Gene Symbol	RGPD5.RGPD6
Uniprot ID	RGPD5_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1560-1640
Immunogen Region	1560-1640
Specificity	RGPD5/RGPD6 polyclonal antibody (Ranbp2-Like And Grip Domain-Containing Protein 5/6) binds to endogenous Ranbp2-Like And Grip Domain-Containing Protein 5/6 at the amino acid region 1560-1640.
Immunogen Sequence	