

Anti-RASGEF1A antibody (130-210) (STJ191262)

STJ191262

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ras-Gef Domain-Containing Family Member 1a (130-210) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	221002
Gene Symbol	RASGEF1A
Uniprot ID	RGF1A_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 130-210
Immunogen Region	130-210
Specificity	RASGEF1A polyclonal antibody (Ras-Gef Domain-Containing Family Member 1a) binds to endogenous Ras-Gef Domain-Containing Family Member 1a at the amino acid region 130-210.
Immunogen Sequence	