

Anti-PTPRU antibody (270-350) (STJ191303)

STJ191303

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Receptor-Type Tyrosine-Protein Phosphatase U (270-350) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10076
Gene Symbol	PTPRU
Uniprot ID	PTPRU_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 270-350
Immunogen Region	270-350
Specificity	PTPRU polyclonal antibody (Receptor-Type Tyrosine-Protein Phosphatase U) binds to endogenous Receptor-Type Tyrosine-Protein Phosphatase U at the amino acid region 270-350.
Immunogen Sequence	