

Anti-ARHGAP6 antibody (480-560) (STJ191324)

STJ191324

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Rho Gtpase-Activating Protein 6 (480-560) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	395
Gene Symbol	ARHGAP6
Uniprot ID	RHG06_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 480-560
Immunogen Region	480-560
Specificity	ARHGAP6 polyclonal antibody (Rho Gtpase-Activating Protein 6) binds to endogenous Rho Gtpase-Activating Protein 6 at the amino acid region 480-560.
Immunogen Sequence	