

Anti-DAGLA antibody (900-980) (STJ191427)

STJ191427

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Diacylglycerol Lipase-Alpha (900-980) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	747
Gene Symbol	DAGLA
Uniprot ID	DGLA_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 900-980
Immunogen Region	900-980
Specificity	DAGLA polyclonal antibody (Diacylglycerol Lipase-Alpha) binds to endogenous Diacylglycerol Lipase-Alpha at the amino acid region 900-980.
Immunogen Sequence	