

Anti-SLC8A2 antibody (380-460) (STJ191438)

STJ191438

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sodium/Calcium Exchanger 2 (380-460) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6543
Gene Symbol	SLC8A2
Uniprot ID	NAC2_HUMAN
Immunogen	Synthesized peptid derived from human protein at aa range 380-460
Immunogen Region	380-460
Specificity	SLC8A2 polyclonal antibody (Sodium/Calcium Exchanger 2) binds to endogenous Sodium/Calcium Exchanger 2 at the amino acid region 380-460.
Immunogen Sequence	