

Anti-TJP1 antibody (STJ191514)

STJ191514

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tight Junction Protein Zo-1 is suitable for use in Western Blot, Immunofluorescence and ELISA research applications.
Applications	WB, IF, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 IF 1:100-300 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7082
Gene Symbol	TJP1
Uniprot ID	ZO1_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Region	
Specificity	TJP1 polyclonal antibody (Tight Junction Protein Zo-1) binds to endogenous Tight Junction Protein Zo-1.
Immunogen Sequence	