

Anti-USP17L1 antibody (STJ191571)

STJ191571

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Ubiquitin Carboxyl-Terminal Hydrolase 17-Like Protein 1 is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB, ELISA |
| Host/Source | Rabbit |
| Reactivity | Human, Rat, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:500-2000 ELISA 1:5000-20000 |
| Formulation | PBS, 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 401447 |
| Gene Symbol | USP17L1 |
| Uniprot ID | U17L1_HUMAN |
| Immunogen | Synthesized peptide derived from part region of human protein |
| Immunogen Region | |
| Specificity | USP17L1 polyclonal antibody (Ubiquitin Carboxyl-Terminal Hydrolase 17-Like Protein 1) binds to endogenous Ubiquitin Carboxyl-Terminal Hydrolase 17-Like Protein 1. |
| Immunogen Sequence | |