

Anti-CDH12 antibody (500-580) (STJ191661)

STJ191661

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cadherin-12 (500-580) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1010
Gene Symbol	CDH12
Uniprot ID	CAD12_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 500-580
Immunogen Region	500-580
Specificity	CDH12 polyclonal antibody (Cadherin-12) binds to endogenous Cadherin-12 at the amino acid region 500-580.
Immunogen Sequence	