

Anti-CCNT2 antibody (510-590) (STJ191670)

STJ191670

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cyclin-T2 (510-590) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	905
Gene Symbol	CCNT2
Uniprot ID	CCNT2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 510-590
Immunogen Region	510-590
Specificity	CCNT2 polyclonal antibody (Cyclin-T2) binds to endogenous Cyclin-T2 at the amino acid region 510-590.
Immunogen Sequence	