

Anti-CUL4A antibody (590-670) (STJ191690)

STJ191690

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Cullin-4a (590-670) is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB, ELISA |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:500-2000 ELISA 1:5000-20000 |
| Formulation | PBS, 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 8451 |
| Gene Symbol | CUL4A |
| Uniprot ID | CUL4A_HUMAN |
| Immunogen | Synthesized peptide derived from human protein at aa range 590-670 |
| Immunogen Region | 590-670 |
| Specificity | CUL4A polyclonal antibody (Cullin-4a) binds to endogenous Cullin-4a at the amino acid region 590-670. |
| Immunogen Sequence | |