

Anti-CUL5 antibody (560-640) (STJ191691)

STJ191691

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cullin-5 (560-640) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8065
Gene Symbol	CUL5
Uniprot ID	CUL5_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 560-640
Immunogen Region	560-640
Specificity	CUL5 polyclonal antibody (Cullin-5) binds to endogenous Cullin-5 at the amino acid region 560-640.
Immunogen Sequence	