

## Anti-MAP4K2 antibody (400-480) (STJ191699)

STJ191699

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Mitogen-Activated Protein Kinase Kinase Kinase Kinase 2 (400-480) is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB, ELISA  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Rat, Mouse  |

### PRODUCT PROPERTIES

|                      |  |
|----------------------|--|
| <b>Clonality</b>     | Polyclonal   |
| <b>Clone ID</b>      |  |
| <b>Concentration</b> | 1 mg/mL  |
| <b>Conjugation</b>   | Unconjugated   |
| <b>Purification</b>  | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.          |
| <b>Dilution</b>      | WB 1:500-2000  |
| <b>Range</b>         | ELISA 1:5000-20000   |
| <b>Formulation</b>   | PBS, 50% Glycerol and 0.02% Sodium Azide.  |
| <b>Isotype</b>       | IgG  |
| <b>Storage</b>       | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
| <b>Instruction</b>   |  |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">5871</a>   |
| <b>Gene Symbol</b>        | <a href="#">MAP4K2</a>   |
| <b>Uniprot ID</b>         | <a href="#">M4K2_HUMAN</a>   |
| <b>Immunogen</b>          | Synthesized peptide derived from human protein at aa range 400-480   |
| <b>Immunogen Region</b>   | 400-480  |
| <b>Specificity</b>        | MAP4K2 polyclonal antibody (Mitogen-Activated Protein Kinase Kinase Kinase Kinase 2) binds to endogenous Mitogen-Activated Protein Kinase Kinase Kinase Kinase 2 at the amino acid region 400-480. |
| <b>Immunogen Sequence</b> |  |