

Anti-KDM4D antibody (400-480) (STJ191746)

STJ191746

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lysine-Specific Demethylase 4d (400-480) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55693
Gene Symbol	KDM4D
Uniprot ID	KDM4D_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 400-480
Immunogen Region	400-480
Specificity	KDM4D polyclonal antibody (Lysine-Specific Demethylase 4d) binds to endogenous Lysine-Specific Demethylase 4d at the amino acid region 400-480.
Immunogen Sequence	