

Anti-LZTR1 antibody (310-390) (STJ191772)

STJ191772

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Leucine-Zipper-Like Transcriptional Regulator 1 (310-390) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8216
Gene Symbol	LZTR1
Uniprot ID	LZTR1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 310-390
Immunogen Region	310-390
Specificity	LZTR1 polyclonal antibody (Leucine-Zipper-Like Transcriptional Regulator 1) binds to endogenous Leucine-Zipper-Like Transcriptional Regulator 1 at the amino acid region 310-390.
Immunogen Sequence	