

Anti-INTS6 antibody (550-630) (STJ191826)

STJ191826

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Integrator Complex Subunit 6 (550-630) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26512
Gene Symbol	INTS6
Uniprot ID	INT6_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 550-630
Immunogen Region	550-630
Specificity	INTS6 polyclonal antibody (Integrator Complex Subunit 6) binds to endogenous Integrator Complex Subunit 6 at the amino acid region 550-630.
Immunogen Sequence	