

Anti-CDY1/CDY1B antibody (210-290) (STJ191832)

STJ191832

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Testis-Specific Chromodomain Protein Y 1 (210-290) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	253175/9085
Gene Symbol	CDY1.CDY1B
Uniprot ID	CDY1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 210-290
Immunogen Region	210-290
Specificity	CDY1/CDY1B polyclonal antibody (Testis-Specific Chromodomain Protein Y 1) binds to endogenous Testis-Specific Chromodomain Protein Y 1 at the amino acid region 210-290.
Immunogen Sequence	