

Anti-FBF1 antibody (730-810) (STJ191845)

STJ191845

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Fas-Binding Factor 1 (730-810) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85302
Gene Symbol	FBF1
Uniprot ID	FBF1_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 730-810
Immunogen Region	730-810
Specificity	FBF1 polyclonal antibody (Fas-Binding Factor 1) binds to endogenous Fas-Binding Factor 1 at the amino acid region 730-810.
Immunogen Sequence	