

Anti-PLCB4 antibody (STJ191892)

STJ191892

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Beta-4 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5332
Gene Symbol	PLCB4
Uniprot ID	PLCB4_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	PLCB4 polyclonal antibody (1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Beta-4) binds to endogenous 1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Beta-4.
Immunogen Sequence	