

Anti-PLCE1 antibody (STJ191893)

STJ191893

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Epsilon-1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51196
Gene Symbol	PLCE1
Uniprot ID	PLCE1_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	PLCE1 polyclonal antibody (1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Epsilon-1) binds to endogenous 1-Phosphatidylinositol 4-5-Bisphosphate Phosphodiesterase Epsilon-1.
Immunogen Sequence	