

Anti-CREB3L3 antibody (STJ191913)

STJ191913

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cyclic Amp-Responsive Element-Binding Protein 3-Like Protein 3 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84699
Gene Symbol	CREB3L3
Uniprot ID	CR3L3_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	CREB3L3 polyclonal antibody (Cyclic Amp-Responsive Element-Binding Protein 3-Like Protein 3) binds to endogenous Cyclic Amp-Responsive Element-Binding Protein 3-Like Protein 3.
Immunogen Sequence	