

Anti-LILRB3 antibody (STJ192018)

STJ192018

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Leukocyte Immunoglobulin-Like Receptor Subfamily B Member 3 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	102725035/107987462/11025
Gene Symbol	LILRB3
Uniprot ID	LIRB3_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	LILRB3 polyclonal antibody (Leukocyte Immunoglobulin-Like Receptor Subfamily B Member 3) binds to endogenous Leukocyte Immunoglobulin-Like Receptor Subfamily B Member 3.
Immunogen Sequence	