

Anti-PPP3CB antibody (STJ192160)

STJ192160

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/Threonine-Protein Phosphatase 2b Catalytic Subunit Beta Isoform is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5532
Gene Symbol	PPP3CB
Uniprot ID	PP2BB_HUMAN
Immunogen	Synthesized peptide derived from part region of human protein
Immunogen Region	
Specificity	PPP3CB polyclonal antibody (Serine/Threonine-Protein Phosphatase 2b Catalytic Subunit Beta Isoform) binds to endogenous Serine/Threonine-Protein Phosphatase 2b Catalytic Subunit Beta Isoform.
Immunogen Sequence	