

Anti-DPP6 antibody (750-830) (STJ192354)

STJ192354

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dipeptidyl Aminopeptidase-Like Protein 6 (750-830) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1804
Gene Symbol	DPP6
Uniprot ID	DPP6_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 750-830
Immunogen Region	750-830
Specificity	DPP6 polyclonal antibody (Dipeptidyl Aminopeptidase-Like Protein 6) binds to endogenous Dipeptidyl Aminopeptidase-Like Protein 6 at the amino acid region 750-830.
Immunogen Sequence	