

Anti-ADAM15 antibody (50-130) (STJ192358)

STJ192358

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Disintegrin And Metalloproteinase Domain-Containing Protein 15 (50-130) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8751
Gene Symbol	ADAM15
Uniprot ID	ADA15_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 50-130
Immunogen Region	50-130
Specificity	ADAM15 polyclonal antibody (Disintegrin And Metalloproteinase Domain-Containing Protein 15) binds to endogenous Disintegrin And Metalloproteinase Domain-Containing Protein 15 at the amino acid region 50-130.
Immunogen Sequence	