

Anti-SEMA3E antibody (570-650) (STJ192382)

STJ192382

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Semaphorin-3e (570-650) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9723
Gene Symbol	SEMA3E
Uniprot ID	SEM3E_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 570-650
Immunogen Region	570-650
Specificity	SEMA3E polyclonal antibody (Semaphorin-3e) binds to endogenous Semaphorin-3e at the amino acid region 570-650.
Immunogen Sequence	