

Anti-PPM1B antibody (310-390) (STJ192389)

STJ192389

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Phosphatase 1b (310-390) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5495
Gene Symbol	PPM1B
Uniprot ID	PPM1B_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 310-390
Immunogen Region	310-390
Specificity	PPM1B polyclonal antibody (Protein Phosphatase 1b) binds to endogenous Protein Phosphatase 1b at the amino acid region 310-390.
Immunogen Sequence	