

Anti-ZSCAN10 antibody (160-240) (STJ192406)

STJ192406

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc Finger And Scan Domain-Containing Protein 10 (160-240) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84891
Gene Symbol	ZSCAN10
Uniprot ID	ZSC10_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 160-240
Immunogen Region	160-240
Specificity	ZSCAN10 polyclonal antibody (Zinc Finger And Scan Domain-Containing Protein 10) binds to endogenous Zinc Finger And Scan Domain-Containing Protein 10 at the amino acid region 160-240.
Immunogen Sequence	