

Anti-MST1 antibody (420-500) (STJ192407)

STJ192407

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hepatocyte Growth Factor-Like Protein (420-500) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4485
Gene Symbol	MST1
Uniprot ID	HGFL_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 420-500
Immunogen Region	420-500
Specificity	MST1 polyclonal antibody (Hepatocyte Growth Factor-Like Protein) binds to endogenous Hepatocyte Growth Factor-Like Protein at the amino acid region 420-500.
Immunogen Sequence	