

Anti-CSPG4 antibody (1910-1990) (STJ192412)

STJ192412

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Chondroitin Sulfate Proteoglycan 4 (1910-1990) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1464
Gene Symbol	CSPG4
Uniprot ID	CSPG4_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 1910-1990
Immunogen Region	1910-1990
Specificity	CSPG4 polyclonal antibody (Chondroitin Sulfate Proteoglycan 4) binds to endogenous Chondroitin Sulfate Proteoglycan 4 at the amino acid region 1910-1990.
Immunogen Sequence	