

Anti-NAALADL2 antibody (190-270) (STJ192421)

STJ192421

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Inactive N-Acetylated-Alpha-Linked Acidic Dipeptidase-Like Protein 2 (190-270) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	254827
Gene Symbol	NAALADL2
Uniprot ID	NADL2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 190-270
Immunogen Region	190-270
Specificity	NAALADL2 polyclonal antibody (Inactive N-Acetylated-Alpha-Linked Acidic Dipeptidase-Like Protein 2) binds to endogenous Inactive N-Acetylated-Alpha-Linked Acidic Dipeptidase-Like Protein 2 at the amino acid region 190-270.
Immunogen Sequence	