

Anti-DISP3 antibody (360-440) (STJ192433)

STJ192433

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein Dispatched Homolog 3 (360-440) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57540
Gene Symbol	DISP3
Uniprot ID	DISP3_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 360-440
Immunogen Region	360-440
Specificity	DISP3 polyclonal antibody (Protein Dispatched Homolog 3) binds to endogenous Protein Dispatched Homolog 3 at the amino acid region 360-440.
Immunogen Sequence	