

Anti-THNSL2 antibody (91-140) (STJ192540)

STJ192540

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Threonine Synthase-Like 2 (91-140) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	PBS, 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55258
Gene Symbol	THNSL2
Uniprot ID	THNS2_HUMAN
Immunogen	Synthesized peptide derived from human protein at aa range 91-140
Immunogen Region	91-140
Specificity	THNSL2 polyclonal antibody (Threonine Synthase-Like 2) binds to endogenous Threonine Synthase-Like 2 at the amino acid region 91-140.
Immunogen Sequence	